

产品名称: USP14 Rabbit Monoclonal Antibody
产品货号: AMRe83952



产品概述 (Summary)

产品名称 (Production Name)	USP14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP14
别名 (Alternative Names)	TGT; tRNA guanine transglycosylase 60 kD subunit; Ubiquitin carboxyl terminal hydrolase 14; Ubiquitin specific peptidase 14; USP14;; USP14
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P54578.A synthesized peptide derived from human USP14

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	56 kDa

研究背景 (Background)


产品名称: USP14 Rabbit Monoclonal Antibody
产品货号: AMRe83952



Ubiquitin-Specific Protease 14, which is also known as the 60 kDa subunit of tRNA-guanine transglycosylase (USP14/TGT60 kDa). USP14 is recruited to the proteasome through its reversible association with the PSMD2 (S2/hRPN1) subunit of the 19S regulatory particle. Whereas PSMD14 appears to promote substrate degradation, USP14 is thought to antagonize substrate degradation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of USP14 expression in (1) HeLa cell lysate; (2) RAW 264.7 cell lysate; (3) C6 cell lysate.

注意事项 (Note)

For research use only .