

产品名称: BNIP1 Rabbit Monoclonal Antibody
产品货号: AMRe83956



产品概述 (Summary)

产品名称 (Production Name)	BNIP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BNIP1 BCL2 adenovirus E1B 19kD interacting protein 1; BNIP1; Nip1; sec20; SEC20L;
别名 (Alternative Names)	Transformation-related gene 8 protein; TRG8; Vesicle transport protein SEC20;;BNIP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q12981.A synthesized peptide derived from human BNIP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	26 kDa

产品名称: BNIP1 Rabbit Monoclonal Antibody
产品货号: AMRe83956




研究背景 (Background)

SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER. Required for maintenance of ER network. Implicated in the suppression of cell death.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of BNIP1 expression in (1) Jurkat cell lysate; (2) C6 cell lysate.

注意事项 (Note)

For research use only .