

产品名称: Sortilin Rabbit Monoclonal Antibody
产品货号: AMRe83958



产品概述 (Summary)

产品名称 (Production Name)	Sortilin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Sortilin
别名 (Alternative Names)	100 kDa NT receptor; Glycoprotein 95; Gp95; LDLCQ6; Neurotensin receptor 3; NT3; NTR3; SORT 1; Sortilin 1;;Sortilin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q99523.A synthesized peptide derived from human Sortilin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	92 kDa

研究背景 (Background)


产品名称: Sortilin Rabbit Monoclonal Antibody
产品货号: AMRe83958




Functions as a sorting receptor in the Golgi compartment and as a clearance receptor on the cell surface. Required for protein transport from the Golgi apparatus to the lysosomes by a pathway that is independent of the mannose-6-phosphate receptor (M6PR). Lysosomal proteins bind specifically to the receptor in the Golgi apparatus and the resulting receptor-ligand complex is transported to an acidic prelysosomal compartment where the low pH mediates the dissociation of the complex.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Sortilin expression in SW480 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: Sortilin Rabbit Monoclonal Antibody
产品货号: AMRe83958



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .