

产品名称: TNF beta Rabbit Monoclonal Antibody
产品货号: AMRe83960



产品概述 (Summary)

产品名称 (Production Name)	TNF beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNF beta
别名 (Alternative Names)	LTalpha; DIF; LT alpha; Lta; Lymphotoxin alpha; TNF B; TNF superfamily member 1; TNF-beta; TNFB; TNFbeta; TNFSF1; TNLG1E;;LT alpha
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01374.A synthesized peptide derived from human LT alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa ; Observed MW: 18 kDa

研究背景 (Background)

产品名称: TNF beta Rabbit Monoclonal Antibody
产品货号: AMRe83960




Cytokine that in its homotrimeric form binds to TNFRSF1A/TNFR1, TNFRSF1B/TNFR and TNFRSF14/HVEM.

In its heterotrimeric form with LTB binds to TNFRSF3/LTBR. Lymphotoxin is produced by lymphocytes and is cytotoxic for a wide range of tumor cells in vitro and in vivo.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TNF beta expression in TNF beta recombinant protein.

注意事项 (Note)

For research use only .