

产品名称: TBPL1 Rabbit Monoclonal Antibody
产品货号: AMRe83964



产品概述 (Summary)

产品名称 (Production Name)	TBPL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBPL1
别名 (Alternative Names)	21 kDa TBP-like protein; Second TBP of unique DNA protein; STUD; TATA box-binding protein-related factor 2; TBPL1; TLF; TLP; TLP21;;TBPL1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P62380.A synthesized peptide derived from human TBPL1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa ; Observed MW: 23 kDa

产品名称: TBPL1 Rabbit Monoclonal Antibody
产品货号: AMRe83964




研究背景 (Background)

Part of a specialized transcription system that mediates the transcription of most ribosomal proteins through the 5'-TCT-3' motif which is a core promoter element at these genes. Seems to also mediate the transcription of NF1. Does not bind the TATA box.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TBPL1 expression in (1) HeLa cell lysate; (1) Mouse testis lysate.

注意事项 (Note)

For research use only .