

产品名称: Golgin97 Rabbit Monoclonal Antibody
产品货号: AMRe83965



产品概述 (Summary)

产品名称 (Production Name)	Golgin97 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Golgin97 Gap junction protein alpha 4 37kD; GOLGA1; Golgi autoantigen golgin
别名 (Alternative Names)	subfamily a 1; Golgin A1; Golgin subfamily A member 1; Golgin97;;Golgin subfamily A member 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92805.A synthesized peptide derived from human Golgin subfamily A member 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 88 kDa ; Observed MW: 97 kDa

产品名称: Golgin97 Rabbit Monoclonal Antibody
产品货号: AMRe83965




研究背景 (Background)

Involved in vesicular trafficking at the Golgi apparatus level. Involved in endosome-to-Golgi trafficking.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Golgin 97 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .