

产品名称: Reptin/RUVBL2 Rabbit Monoclonal Antibody
产品货号: AMRe83967



产品概述 (Summary)

产品名称 (Production Name)	Reptin/RUVBL2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Reptin/RUVBL2
别名 (Alternative Names)	ECP51; Reptin 52; REPTIN; RUVB2; RUVBL2; RVB2; TAP54-beta; TIH2; TIP48; TIP49b; zreptin;;RUVBL2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y230.A synthesized peptide derived from human RUVBL2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 51 kDa ; Observed MW: 48 kDa

研究背景 (Background)


产品名称: Reptin/RUVBL2 Rabbit Monoclonal Antibody
产品货号: AMRe83967



Possesses single-stranded DNA-stimulated ATPase and ATP-dependent DNA helicase (5' to 3') activity; hexamerization is thought to be critical for ATP hydrolysis and adjacent subunits in the ring-like structure contribute to the ATPase activity. Component of the NuA4 histone acetyltransferase complex which is involved in transcriptional activation of select genes principally by acetylation of nucleosomal histones H4 and H2A.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Reptin / RUVBL2 expression in (1) HeLa cell lysate; (2) NIH/3T3 cell lysate.

注意事项 (Note)

For research use only .