

产品名称: Prothrombin Rabbit Monoclonal Antibody
产品货号: AMRe83968



产品概述 (Summary)

产品名称 (Production Name)	Prothrombin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Prothrombin Coagulation factor II; F2; Factor II; Prepro coagulation factor II; Prothrombin;
别名 (Alternative Names)	prothrombin B-chain; PT; RPRGL2; serine protease; THPH1; Thrombin heavy chain;;Prothrombin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P00734.A synthesized peptide derived from human Prothrombin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	70 kDa

产品名称: Prothrombin Rabbit Monoclonal Antibody
产品货号: AMRe83968




研究背景 (Background)

Thrombin, which cleaves bonds after Arg and Lys, converts fibrinogen to fibrin and activates factors V, VII, VIII, XIII, and, in complex with thrombomodulin, protein C. Functions in blood homeostasis, inflammation and wound healing.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Prothrombin expression in human serum lysate.

注意事项 (Note)

For research use only .