

产品名称: ITCH Rabbit Monoclonal Antibody
产品货号: AMRe83971



产品概述 (Summary)

产品名称 (Production Name)	ITCH Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITCH ADMFD; AIF4; AIP4; Atrophin 1 interacting protein 4; Itchy homolog E3
别名 (Alternative Names)	ubiquitin protein ligase; NAPP1; NFE2 associated polypeptide 1; Ubiquitin protein ligase ITCH;;ITCH
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96J02.A synthesized peptide derived from human ITCH

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	103 kDa

产品名称: ITCH Rabbit Monoclonal Antibody
产品货号: AMRe83971




研究背景 (Background)


Acts as an E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of ITCH expression in K562 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: ITCH Rabbit Monoclonal Antibody
产品货号: AMRe83971



注意事项 (Note)

For research use only .