

产品名称: ORM1 Rabbit Monoclonal Antibody
产品货号: AMRe83975



产品概述 (Summary)

产品名称 (Production Name)	ORM1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ORM1
别名 (Alternative Names)	AGP 1; AGP A; AGP1; Alpha 1 acid glycoprotein; alpha-1-AGP; HEL S 153w; ORM; Orosomucoid 1;;AGP 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P02763.A synthesized peptide derived from human AGP 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa ; Observed MW: 47 kDa

研究背景 (Background)


产品名称: ORM1 Rabbit Monoclonal Antibody
产品货号: AMRe83975



Functions as a transport protein in the blood stream. Binds various ligands in the interior of its beta-barrel domain. Also binds synthetic drugs and influences their distribution and availability in the body. Appears to function in modulating the activity of the immune system during the acute-phase reaction.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ORM1 expression in human testis lysate.

注意事项 (Note)

For research use only .