

产品名称: INDO Rabbit Monoclonal Antibody
产品货号: AMRe83976



产品概述 (Summary)

产品名称 (Production Name)	INDO Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	INDO
别名 (Alternative Names)	3-dioxygenase; IDO 1; IDO; IDO1; INDO; indolamine 2,3 dioxygenase; Indole 2 3 dioxygenase;; IDO 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P14902.A synthesized peptide derived from human IDO 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	45 kDa

研究背景 (Background)


产品名称: INDO Rabbit Monoclonal Antibody
产品货号: AMRe83976



Catalyzes the first and rate limiting step of the catabolism of the essential amino acid tryptophan along the kynurenine pathway. Involved in the peripheral immune tolerance, contributing to maintain homeostasis by preventing autoimmunity or immunopathology that would result from uncontrolled and overreacting immune responses.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of INDO expression in HeLa cell lysate treated with IFN gamma.

注意事项 (Note)

For research use only .