

产品名称: KIF3A Rabbit Monoclonal Antibody
产品货号: AMRe83977



产品概述 (Summary)

产品名称 (Production Name)	KIF3A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KIF3A KIF 3; KIF 3A; Kifl; Kinesin family member 3A; Kinesin like protein KIF 3A;
别名 (Alternative Names)	Kinesin-like protein KIF3A; Kns3; ;KIF3A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y496.A synthesized peptide derived from human KIF3A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	80 kDa

产品名称: KIF3A Rabbit Monoclonal Antibody
产品货号: AMRe83977




研究背景 (Background)

Microtubule-based anterograde translocator for membranous organelles. Plus end-directed microtubule sliding activity in vitro. Plays a role in primary cilia formation. Plays a role in centriole cohesion and subdistal appendage organization and function. Regulates the formation of the subdistal appendage via recruitment of DCTN1 to the centriole. Also required for ciliary basal feet formation and microtubule anchoring to mother centriole.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of KIF3A expression in SH-SY5Y cell lysate.

注意事项 (Note)

For research use only .