

产品名称: GLUR Rabbit Monoclonal Antibody
产品货号: AMRe83978



产品概述 (Summary)

产品名称 (Production Name)	GLUR Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GLUR GLR1; GLUR; Glutathione reductase; Glutathione reductase mitochondrial;
别名 (Alternative Names)	Glutathione reductase, mitochondrial; GR; Gr1; GRase; GRD1; GSR; ;Glutathione reductase
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P00390.A synthesized peptide derived from human Glutathione reductase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	56 kDa

产品名称: GLUR Rabbit Monoclonal Antibody
产品货号: AMRe83978




研究背景 (Background)


Maintains high levels of reduced glutathione in the cytosol.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GLUR expression in (1) HeLa cell lysate; (2) Mouse brain lysate; (3) Rat kidney lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .