

产品名称: DDX58 Rabbit Monoclonal Antibody
产品货号: AMRe83979



产品概述 (Summary)

产品名称 (Production Name)	DDX58 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDX58
别名 (Alternative Names)	RIG-I-like receptor 1; RLR-1; RLR1; Retinoic acid-inducible gene 1 protein; RIG-1; RIG1; RIG-I; RIGI; RIGI
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95786.A synthesized peptide derived from human RIGI

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	107 kDa

研究背景 (Background)


产品名称: DDX58 Rabbit Monoclonal Antibody
产品货号: AMRe83979



Innate immune receptor that senses cytoplasmic viral nucleic acids and activates a downstream signaling cascade leading to the production of type I interferons and pro-inflammatory cytokines.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DDX58 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .