

产品名称: Grp75 Rabbit Monoclonal Antibody
产品货号: AMRe83981



产品概述 (Summary)

产品名称 (Production Name)	Grp75 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Grp75
别名 (Alternative Names)	GRP-75; Heat shock 70 kDa protein 9; Mortalin; MOT; Peptide-binding protein 74; PBP74; HSPA9; GRP75; HSPA9B; mt-HSP70;;Hsp70 9
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P38646.A synthesized peptide derived from human Hsp70 9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	74 kDa

研究背景 (Background)


产品名称: Grp75 Rabbit Monoclonal Antibody
产品货号: AMRe83981




Chaperone protein which plays an important role in mitochondrial iron-sulfur cluster (ISC) biogenesis. Interacts with and stabilizes ISC cluster assembly proteins FXN, NFU1, NFS1 and ISCU. Regulates erythropoiesis via stabilization of ISC assembly.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Grp75 expression in (1) HeLa cell lysate; (2) Mouse liver lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: Grp75 Rabbit Monoclonal Antibody
产品货号: AMRe83981



注意事项 (Note)

For research use only .