

产品名称: Nuclear Matrix Protein p84 Rabbit Monoclonal Antibody
产品货号: AMRe83984



产品概述 (Summary)

产品名称 (Production Name)	Nuclear Matrix Protein p84 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nuclear Matrix Protein p84
别名 (Alternative Names)	Death domain containing protein p84N5; HPR1; hTRES84; Nuclear matrix protein p84; p84N5; THO complex 1; Tho1; Thoc1;;THOC1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96FV9.A synthesized peptide derived from human THOC1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa ; Observed MW: 84 kDa

产品名称: Nuclear Matrix Protein p84 Rabbit Monoclonal Antibody
产品货号: AMRe83984




研究背景 (Background)

Required for efficient export of polyadenylated RNA. Acts as component of the THO subcomplex of the TREX complex which is thought to couple mRNA transcription, processing and nuclear export, and which specifically associates with spliced mRNA and not with unspliced pre-mRNA.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Nuclear Matrix Protein p84 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .