

**产品名称: ProteinC Rabbit Monoclonal Antibody**  
**产品货号: AMRe83988**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | ProteinC Rabbit Monoclonal Antibody    |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB                                     |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                   |
| 形式 (Form)           | Liquid                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.     |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                                          |
|------------------------|------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | ProteinC<br>Activation peptide; Anticoagulant protein C; APC; Autoprothrombin IIA; Blood |
| 别名 (Alternative Names) | coagulation factor XIV; PROC1; THPH3; THPH4; Vitamin K dependent protein C; ;PROC        |
| 基因 ID (Gene ID)        |                                                                                          |
| 蛋白 ID (SwissProt ID)   | P04070.A synthesized peptide derived from human PROC                                     |

## 产品应用 (Application)

|                          |                  |
|--------------------------|------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000 |
| 蛋白分子量 (Molecular Weight) | 52 kDa           |

**产品名称: ProteinC Rabbit Monoclonal Antibody**  
**产品货号: AMRe83988**




## 研究背景 (Background)

Protein C is a vitamin K-dependent serine protease that regulates blood coagulation by inactivating factors Va and VIIIa in the presence of calcium ions and phospholipids. Exerts a protective effect on the endothelial cell barrier function.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of Protein C expression in HepG2 cell lysate.

## 注意事项 (Note)

For research use only .