

产品名称: Dynein intermediate chain 1 Rabbit Monoclonal Antibody  
产品货号: AMRe83991

---

## 产品概述 (Summary)

|                        |                                                        |
|------------------------|--------------------------------------------------------|
| 产品名称 (Production Name) | Dynein intermediate chain 1 Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody                 |
| 宿主 (Host)              | Rabbit                                                 |
| 应用 (Application)       | WB                                                     |
| 种属反应性 (Reactivity)     | Human, Mouse, Rat                                      |

## 产品性能 (Performance)

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                   |
| 形式 (Form)           | Liquid                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.     |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                |
|------------------------|----------------------------------------------------------------|
| 基因名 (Gene Name)        | Dynein intermediate chain 1                                    |
| 别名 (Alternative Names) | DNAI1; ICS; Immotile cilia syndrome 1; ICS1; CILD1; PCD;;DNAI1 |
| 基因 ID (Gene ID)        |                                                                |
| 蛋白 ID (SwissProt ID)   | Q9UI46.A synthesized peptide derived from human DNAI1          |

## 产品应用 (Application)

|                          |                  |
|--------------------------|------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000 |
| 蛋白分子量 (Molecular Weight) | 79 kDa           |

## 研究背景 (Background)

---


**产品名称: Dynein intermediate chain 1 Rabbit Monoclonal Antibody**  
**产品货号: AMRe83991**



Part of the dynein complex of respiratory cilia.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of Dynein intermediate chain 1 expression in mouse testis lysate.

## 注意事项 (Note)

For research use only .