

产品名称: Superoxide Dismutase 3 Rabbit Monoclonal Antibody
产品货号: AMRe83992



产品概述 (Summary)

产品名称 (Production Name)	Superoxide Dismutase 3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Superoxide Dismutase 3
别名 (Alternative Names)	EC SOD; Extracellular superoxide dismutase [Cu Zn]; Extracellular superoxide dismutase precursor; SOD3; Superoxide dismutase 3 extracellular; ;SOD3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P08294.A synthesized peptide derived from human SOD3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	26 kDa

产品名称: Superoxide Dismutase 3 Rabbit Monoclonal Antibody
产品货号: AMRe83992




研究背景 (Background)


Protect the extracellular space from toxic effect of reactive oxygen intermediates by converting superoxide radicals into hydrogen peroxide and oxygen.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Superoxide Dismutase 3 expression in human fetal heart lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .