

产品名称: MUC4 Rabbit Monoclonal Antibody
产品货号: AMRe83995



产品概述 (Summary)

产品名称 (Production Name)	MUC4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MUC4 Ascites sialoglycoprotein 1; Ascites sialoglycoprotein 2; Ascites sialoglycoprotein; ASGP; ASGP-1; ASGP-2; Muc4; Mucin 4 cell surface associated; Mucin 4 tracheobronchial; Mucin-4 beta chain; Testis mucin;; Mucin 4
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q99102.A synthesized peptide derived from human Mucin 4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 542 kDa ; Observed MW: 120 kDa

产品名称: MUC4 Rabbit Monoclonal Antibody
产品货号: AMRe83995




研究背景 (Background)


Membrane-bound mucin, a family of highly glycosylated proteins that constitute the major component of the mucus, the slimy and viscous secretion covering epithelial surfaces. These glycoproteins play important roles in the protection of the epithelium and are implicated in epithelial renewal and differentiation.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of MUC4 expression in MCF7 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: MUC4 Rabbit Monoclonal Antibody
产品货号: AMRe83995



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .