

产品名称: USP28 Rabbit Monoclonal Antibody
产品货号: AMRe84008



产品概述 (Summary)

产品名称 (Production Name)	USP28 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP28
别名 (Alternative Names)	KIAA1515; Ubiquitin carboxyl terminal hydrolase 28; Ubiquitin thioesterase 28; UBP28;; USP28
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96RU2.A synthesized peptide derived from human USP28

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 122 kDa ; Observed MW: 135 kDa

研究背景 (Background)


产品名称: USP28 Rabbit Monoclonal Antibody
产品货号: AMRe84008



Deubiquitinase involved in DNA damage response checkpoint and MYC proto-oncogene stability. Involved in DNA damage induced apoptosis by specifically deubiquitinating proteins of the DNA damage pathway such as CLSPN. Also involved in G2 DNA damage checkpoint, by deubiquitinating CLSPN, and preventing its degradation by the anaphase promoting complex/cyclosome (APC/C).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of USP28 expression in A431 cell lysate.

注意事项 (Note)

For research use only .