

产品名称: RNF8 Rabbit Monoclonal Antibody
产品货号: AMRe84009



产品概述 (Summary)

产品名称 (Production Name)	RNF8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF8
别名 (Alternative Names)	C3HC4 type zinc finger protein; RNF8;;RNF8
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O76064.A synthesized peptide derived from human RNF8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa ; Observed MW: 60 kDa

研究背景 (Background)

E3 ubiquitin-protein ligase that plays a key role in DNA damage signaling via 2 distinct roles: by mediating the 'Lys-63'-


产品名称: RNF8 Rabbit Monoclonal Antibody
产品货号: AMRe84009



linked ubiquitination of histones H2A and H2AX and promoting the recruitment of DNA repair proteins at double-strand breaks (DSBs) sites, and by catalyzing 'Lys-48'-linked ubiquitination to remove target proteins from DNA damage sites.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RNF8 expression in human fetal kidney lysate.

注意事项 (Note)

For research use only .