

产品名称: TAX1BP3 Rabbit Monoclonal Antibody
产品货号: AMRe84010



产品概述 (Summary)

产品名称 (Production Name)	TAX1BP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAX1BP3
别名 (Alternative Names)	Glutaminase interacting protein 3; Tax interaction protein 1; Tax1 binding protein 3; Tax1bp3; TIP1;;Tax1 binding protein 3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O14907.A synthesized peptide derived from human Tax1 binding protein 3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	14 kDa

研究背景 (Background)


产品名称: TAX1BP3 Rabbit Monoclonal Antibody
产品货号: AMRe84010



May regulate a number of protein-protein interactions by competing for PDZ domain binding sites. Binds CTNNB1 and may thereby act as an inhibitor of the Wnt signaling pathway. Competes with LIN7A for KCNJ4 binding, and thereby promotes KCNJ4 internalization. May play a role in the Rho signaling pathway. May play a role in activation of CDC42 by the viral protein HPV16 E6.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TAX1BP3 expression in HT-1080 cell lysate.

注意事项 (Note)

For research use only .