

产品名称: PI3 Kinase p110 delta Rabbit Monoclonal Antibody
产品货号: AMRe84012



产品概述 (Summary)

产品名称 (Production Name)	PI3 Kinase p110 delta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PI3 Kinase p110 delta
别名 (Alternative Names)	APDS; GRB1; IMD14; p110dp85a; p85-ALPHA; Phosphoinositide 3 kinase B; Phosphoinositide 3 kinase C; p37delta; PI3 kinase p110 subunit delta; PI3Kdelta; Pik3cd; PIK3R1; PK3CD; PtdIns 3 kinase p110;; PI3 kinase p110 delta
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O00329.A synthesized peptide derived from human PI3 kinase p110 delta

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	119 kDa

产品名称: PI3 Kinase p110 delta Rabbit Monoclonal Antibody
产品货号: AMRe84012




研究背景 (Background)


Phosphoinositide-3-kinase (PI3K) phosphorylates phosphatidylinositol (PI) and its phosphorylated derivatives at position 3 of the inositol ring to produce 3-phosphoinositides.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.


产品名称: PI3 Kinase p110 delta Rabbit Monoclonal Antibody
产品货号: AMRe84012



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

产品名称: PI3 Kinase p110 delta Rabbit Monoclonal Antibody
产品货号: AMRe84012



Immunofluorescent analysis using the Antibody at 1:50 dilution.

注意事项 (Note)

For research use only .