

产品名称: NuMA Rabbit Monoclonal Antibody
产品货号: AMRe84025



产品概述 (Summary)

产品名称 (Production Name)	NuMA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NuMA
别名 (Alternative Names)	NMP 22; NUMA; NUMA1; SP H antigen;;NUMA1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q14980.A synthesized peptide derived from human NUMA1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	238 kDa

研究背景 (Background)

Microtubule (MT)-binding protein that plays a role in the formation and maintenance of the spindle poles and the


产品名称: NuMA Rabbit Monoclonal Antibody
产品货号: AMRe84025




alignment and the segregation of chromosomes during mitotic cell division. Functions to tether the minus ends of MTs at the spindle poles, which is critical for the establishment and maintenance of the spindle poles

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Mad2L1 expression in HeLa cell lysate.



All lanes use the Antibody at 1:4K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:4K dilution for 1 hour at room temperature.

产品名称: NuMA Rabbit Monoclonal Antibody
产品货号: AMRe84025



注意事项 (Note)

For research use only .