

产品名称: PFDN5 Rabbit Monoclonal Antibody
产品货号: AMRe84027



产品概述 (Summary)

产品名称 (Production Name)	PFDN5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PFDN5
别名 (Alternative Names)	EIG 1; Eig1; MM1; PFD5; PFDN5; Prefoldin 5; ;PFDN5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q99471.A synthesized peptide derived from human PFDN5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	17 kDa

研究背景 (Background)

Binds specifically to cytosolic chaperonin (c-CPN) and transfers target proteins to it. Binds to nascent polypeptide chain and


产品名称: PFDN5 Rabbit Monoclonal Antibody
产品货号: AMRe84027



promotes folding in an environment in which there are many competing pathways for nonnative proteins. Represses the transcriptional activity of MYC.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PFDN5 expression in 293T cell lysate.

注意事项 (Note)

For research use only .