

产品名称: WIPI1 Rabbit Monoclonal Antibody
产品货号: AMRe84029



产品概述 (Summary)

产品名称 (Production Name)	WIPI1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WIPI1
别名 (Alternative Names)	ATG18; ATG18A; WIPI 1; WIPI-1; WIPI49;;WIPI 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q5MNZ9.A synthesized peptide derived from human WIPI 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 49 kDa ; Observed MW: 48 kDa

研究背景 (Background)

Component of the autophagy machinery that controls the major intracellular degradation process by which cytoplasmic


产品名称: WIPI1 Rabbit Monoclonal Antibody
产品货号: AMRe84029




materials are packaged into autophagosomes and delivered to lysosomes for degradation.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of WIPI1 expression in A375 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: WIPI1 Rabbit Monoclonal Antibody
产品货号: AMRe84029



注意事项 (Note)

For research use only .