

产品名称: ATP5G1 Rabbit Monoclonal Antibody
产品货号: AMRe84035



产品概述 (Summary)

产品名称 (Production Name)	ATP5G1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATP5G1
别名 (Alternative Names)	ATP5G1; ATP5G2; ATP5G3;;ATP5MC1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P05496.A synthesized peptide derived from human ATP5MC1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 14 kDa ; Observed MW: 8 kDa

研究背景 (Background)

Mitochondrial membrane ATP synthase (F(1)F(0) ATP synthase or Complex V) produces ATP from ADP in the presence of a


产品名称: ATP5G1 Rabbit Monoclonal Antibody
产品货号: AMRe84035



proton gradient across the membrane which is generated by electron transport complexes of the respiratory chain.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ATP5G1 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .