

产品名称: CECR5 Rabbit Monoclonal Antibody
产品货号: AMRe84045



产品概述 (Summary)

产品名称 (Production Name)	CECR5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CECR5
别名 (Alternative Names)	Cecr5;;HDHD5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BXW7.A synthesized peptide derived from human HDHD5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	46 kDa

研究背景 (Background)

Candidate gene for the Cat Eye Syndrome (CES), a developmental disorder associated with the duplication of a 2 Mb region


产品名称: CECR5 Rabbit Monoclonal Antibody
产品货号: AMRe84045




of 22q11.2. Duplication usually takes in the form of a supernumerary bisatellited isodicentric chromosome, resulting in four copies of the region (represents an inv dup22(q11)). CES is characterized clinically by the combination of coloboma of the iris and anal atresia with fistula, downslanting palpebral fissures, preauricular tags and/or pits, frequent occurrence of heart and renal malformations, and normal or near-normal mental development.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CECR5 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1k dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .