

产品名称: MSRB3 Rabbit Monoclonal Antibody
产品货号: AMRe84049



产品概述 (Summary)

产品名称 (Production Name)	MSRB3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MSRB3
别名 (Alternative Names)	DFNB74; MsrB3;;MSRB3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8IXL7.A synthesized peptide derived from human MSRB3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa ; Observed MW: 18 kDa

研究背景 (Background)


Catalyzes the reduction of free and protein-bound methionine sulfoxide to methionine. Isoform 2 is essential for hearing.

产品名称: MSRB3 Rabbit Monoclonal Antibody
产品货号: AMRe84049




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of MSRB3 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

产品名称: MSRB3 Rabbit Monoclonal Antibody
产品货号: AMRe84049



For research use only .