

产品名称: PTPRR Rabbit Monoclonal Antibody
产品货号: AMRe84050



产品概述 (Summary)

产品名称 (Production Name)	PTPRR Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PTPRR
别名 (Alternative Names)	ECPTP; NC PTPCOM1; PTPBR7; PTPRQ; PTPRR; PTPSL; Receptor type tyrosine protein phosphatase R;;PTPRR
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15256.A synthesized peptide derived from human PTPRR

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	74 kDa

研究背景 (Background)


产品名称: PTPRR Rabbit Monoclonal Antibody
产品货号: AMRe84050




Sequesters mitogen-activated protein kinases (MAPKs) such as MAPK1, MAPK3 and MAPK14 in the cytoplasm in an inactive form. The MAPKs bind to a dephosphorylated kinase interacting motif, phosphorylation of which by the protein kinase A complex releases the MAPKs for activation and translocation into the nucleus (By similarity).

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of PTPRR expression in mouse brain lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: PTPRR Rabbit Monoclonal Antibody
产品货号: AMRe84050



注意事项 (Note)

For research use only .