

产品名称: SAFB Rabbit Monoclonal Antibody
产品货号: AMRe84056



产品概述 (Summary)

产品名称 (Production Name)	SAFB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SAFB
别名 (Alternative Names)	HAP; HET; SAFB; SAFB1;;SAF B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15424.A synthesized peptide derived from human SAF B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 103 kDa ; Observed MW: 150 kDa

研究背景 (Background)

Binds to scaffold/matrix attachment region (S/MAR) DNA and forms a molecular assembly point to allow the formation of a


产品名称: SAFB Rabbit Monoclonal Antibody
产品货号: AMRe84056



'transcriptosomal' complex (consisting of SR proteins and RNA polymerase II) coupling transcription and RNA processing. Functions as an estrogen receptor corepressor and can also bind to the HSP27 promoter and decrease its transcription.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SAFB expression in 293 cell lysate.

注意事项 (Note)

For research use only .