

产品名称: AZI2 Rabbit Monoclonal Antibody
产品货号: AMRe84059



产品概述 (Summary)

产品名称 (Production Name)	AZI2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AZI2
别名 (Alternative Names)	AZ2; AZI2; Nak associated protein 1; NAP1; TILP;;AZI2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9H6S1.A synthesized peptide derived from human AZI2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	45 kDa

研究背景 (Background)

Adapter protein which binds TBK1 and IKKε playing a role in antiviral innate immunity. Activates serine/threonine-protein


产品名称: AZI2 Rabbit Monoclonal Antibody
产品货号: AMRe84059




kinase TBK1 and facilitates its oligomerization. Enhances the phosphorylation of NF-kappa-B p65 subunit RELA by TBK1.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of AZI2 expression in HeLa cell lysate.

产品名称: AZI2 Rabbit Monoclonal Antibody
产品货号: AMRe84059



Immunofluorescent analysis using the Antibody at 1:50 dilution.

注意事项 (Note)

For research use only .