

产品名称: DNAJA4 Rabbit Monoclonal Antibody
产品货号: AMRe84074



产品概述 (Summary)

产品名称 (Production Name)	DNAJA4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAJA4
别名 (Alternative Names)	DNAJA4; MST104; MSTP104; PRO1472;;DNAJA4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8WW22.A synthesized peptide derived from human DNAJA4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	45 kDa

研究背景 (Background)

Enables chaperone binding activity and unfolded protein binding activity. Involved in several processes, including negative


产品名称: DNAJA4 Rabbit Monoclonal Antibody
产品货号: AMRe84074




regulation of endothelial cell migration; negative regulation of inclusion body assembly; and protein refolding. Located in cytosol and membrane.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of DNAJA4 expression in HeLa cell lysate.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: DNAJA4 Rabbit Monoclonal Antibody
产品货号: AMRe84074



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:500 dilution)

注意事项 (Note)

For research use only .