

产品名称: PSMB10 Rabbit Monoclonal Antibody
产品货号: AMRe84078



产品概述 (Summary)

产品名称 (Production Name)	PSMB10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PSMB10
别名 (Alternative Names)	LMP10; MECL1; Psmb10;; PSMB10
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P40306.A synthesized peptide derived from human PSMB10

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa ; Observed MW: 27 kDa

研究背景 (Background)

The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe,


产品名称: PSMB10 Rabbit Monoclonal Antibody
产品货号: AMRe84078



Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. The proteasome has an ATP-dependent proteolytic activity. This subunit is involved in antigen processing to generate class I binding peptides.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PSMB10 expression in Daudi cell lysate.

注意事项 (Note)

For research use only .