

产品名称: NDUFS2 Rabbit Monoclonal Antibody
产品货号: AMRe84080



产品概述 (Summary)

产品名称 (Production Name)	NDUFS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NDUFS2
别名 (Alternative Names)	Ndufs2;;NDUFS2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75306.A synthesized peptide derived from human NDUFS2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa ; Observed MW: 45 kDa

研究背景 (Background)

Core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I) which catalyzes electron


产品名称: NDUFS2 Rabbit Monoclonal Antibody
产品货号: AMRe84080



transfer from NADH through the respiratory chain, using ubiquinone as an electron acceptor.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NDUFS2 expression in (1) HeLa cell lysate; (2) RAW264.7 HeLa cell lysate.

注意事项 (Note)

For research use only .