

产品名称: UCHL5IP Rabbit Monoclonal Antibody
产品货号: AMRe84081



产品概述 (Summary)

产品名称 (Production Name)	UCHL5IP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UCHL5IP
别名 (Alternative Names)	Haus7; HSXQ28ORF; UCHL5IP; UIP1;;HAUS7
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q99871.A synthesized peptide derived from human HAUS7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa ; Observed MW: 37 kDa

研究背景 (Background)

Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the


产品名称: UCHL5IP Rabbit Monoclonal Antibody
产品货号: AMRe84081




HAUS augmin-like complex.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of UCHL5IP expression in Jurkat cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: UCHL5IP Rabbit Monoclonal Antibody
产品货号: AMRe84081



注意事项 (Note)

For research use only .