

产品名称: PRPF4 Rabbit Monoclonal Antibody
产品货号: AMRe84088



产品概述 (Summary)

产品名称 (Production Name)	PRPF4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRPF4
别名 (Alternative Names)	hPrp4; HPRP4P; PRP4; Prp4p; PRPF4; RP70; SNRNP60;;PRPF4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43172.A synthesized peptide derived from human PRPF4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	58 kDa

研究背景 (Background)

Plays a role in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex that is involved in spliceosome


产品名称: PRPF4 Rabbit Monoclonal Antibody
产品货号: AMRe84088



assembly, and as component of the precatalytic spliceosome (spliceosome B complex).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PRPF4 expression in (1) Raji cell lysate; (2) RAW 264.7 cell lysate; (3) C6 cell lysate.

注意事项 (Note)

For research use only .