

产品名称: ATE1 Rabbit Monoclonal Antibody
产品货号: AMRe84089



产品概述 (Summary)

产品名称 (Production Name)	ATE1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATE1
别名 (Alternative Names)	Arginyltransferase 1; Ate1; R transferase 1;;Arginyltransferase 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O95260.A synthesized peptide derived from human Arginyltransferase 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	59 kDa

研究背景 (Background)

Involved in the post-translational conjugation of arginine to the N-terminal aspartate or glutamate of a protein. This


产品名称: ATE1 Rabbit Monoclonal Antibody
产品货号: AMRe84089



arginylation is required for degradation of the protein via the ubiquitin pathway. Does not arginylate cysteine residues.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ATE1 expression in (1) HepG2 cell lysate; (2) Mouse spleen lysate.

注意事项 (Note)

For research use only .