

产品名称: PMF1 Rabbit Monoclonal Antibody
产品货号: AMRe84090



产品概述 (Summary)

产品名称 (Production Name)	PMF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PMF1
别名 (Alternative Names)	Pmf1; Polyamine modulated factor 1;;Polyamine modulated factor 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q6P1K2.A synthesized peptide derived from human Polyamine modulated factor 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa ; Observed MW: 25 kDa

研究背景 (Background)


产品名称: PMF1 Rabbit Monoclonal Antibody
产品货号: AMRe84090



Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochore formation during mitosis. May act as a cotranscription partner of NFE2L2 involved in regulation of polyamine-induced transcription of SSAT.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PMF1 expression in (1) HeLa cell lysate; (2) RAW 264.7 cell lysate.

注意事项 (Note)

For research use only .