

产品名称: Adiponectin Rabbit Monoclonal Antibody
产品货号: AMRe84091



产品概述 (Summary)

产品名称 (Production Name)	Adiponectin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Adiponectin 30 kDa adipocyte complement related protein; ACDC; ACRP30; Adipocyte; Adipocyte C1q and collagen domain containing protein; Adipocyte complement related 30 kDa protein; Adipocyte complement related protein of 30 kDa; adipocyte-specific secretory protein; Adiponectin; Adiponectin precursor; adiponectin, C1Q and collagen domain containing; Adipoq; Adipose most abundant gene transcript 1; Adipose specific collagen like factor; ADIPQTL1; ADPN; APM1; C1q and collagen domain-containing protein; GBP28; Gelatin binding protein; Gelatin binding protein 28; ;Adiponectin
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q60994(mouse)(human:Q15848).A synthesized peptide derived from mouse Adiponectin

产品名称: Adiponectin Rabbit Monoclonal Antibody
产品货号: AMRe84091



产品应用 (Application)


稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200, ICC/IF 1:50-1:200, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 27 kDa ; Observed MW: 28 kDa

研究背景 (Background)

Important adipokine involved in the control of fat metabolism and insulin sensitivity, with direct anti-diabetic, anti-atherogenic and anti-inflammatory activities. Stimulates AMPK phosphorylation and activation in the liver and the skeletal muscle, enhancing glucose utilization and fatty-acid combustion.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Adiponectin expression in mouse kidney cell lysate.

注意事项 (Note)

For research use only .