

产品名称: SSX2IP Rabbit Monoclonal Antibody
产品货号: AMRe84092



产品概述 (Summary)

产品名称 (Production Name)	SSX2IP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SSX2IP
别名 (Alternative Names)	ADIP; SSX2IP;;SSX2IP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y2D8.A synthesized peptide derived from human SSX2IP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 71 kDa ; Observed MW: 74,76 kDa

研究背景 (Background)

Belongs to an adhesion system, which plays a role in the organization of homotypic, interneuronal and heterotypic cell-cell


产品名称: SSX2IP Rabbit Monoclonal Antibody
产品货号: AMRe84092



adherens junctions (AJs). May connect the nectin-afadin and E-cadherin-catenin system through alpha-actinin and may be involved in organization of the actin cytoskeleton at AJs through afadin and alpha-actinin. Involved in cell movement: localizes at the leading edge of moving cells in response to PDGF and is required for the formation of the leading edge and the promotion of cell movement, possibly via activation of Rac signaling.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SSX2IP expression in HeLa cell lysate.

注意事项 (Note)

For research use only .