

产品名称: SLC25A12 Rabbit Monoclonal Antibody
产品货号: AMRe84093



产品概述 (Summary)

产品名称 (Production Name)	SLC25A12 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC25A12
别名 (Alternative Names)	AGC1; ARALAR; ARALAR1; SLC25A12;;SLC25A12
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75746.A synthesized peptide derived from human SLC25A12

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 75 kDa ; Observed MW: 72 kDa

研究背景 (Background)

Mitochondrial electrogenic aspartate/glutamate antiporter that favors efflux of aspartate and entry of glutamate and


产品名称: SLC25A12 Rabbit Monoclonal Antibody
产品货号: AMRe84093



proton within the mitochondria as part of the malate-aspartate shuttle.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SLC25A12 expression in (1) A431 cell lysate; (2) NIH/3T3 cell lysate.

注意事项 (Note)

For research use only .