

产品名称: RPS8 Rabbit Monoclonal Antibody
产品货号: AMRe84096



产品概述 (Summary)

产品名称 (Production Name)	RPS8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS8
别名 (Alternative Names)	RPS8;;RPS8
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P62241.A synthesized peptide derived from human RPS8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	24 kDa

研究背景 (Background)


Belongs to the ribosomal protein S8e family.

产品名称: RPS8 Rabbit Monoclonal Antibody
产品货号: AMRe84096




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RPS8 expression in (1) 293 cell lysate; (2) RAW 264.7 cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:1K dilution)

注意事项 (Note)

For research use only .