

产品名称: Annexin A7 Rabbit Monoclonal Antibody
产品货号: AMRe84099



产品概述 (Summary)

产品名称 (Production Name)	Annexin A7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Annexin A7
别名 (Alternative Names)	Annexin A7; Annexin VII; ANX7; ANXA7; SNX; Synexin;;Annexin A7
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P20073.A synthesized peptide derived from human Annexin A7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa ; Observed MW: 45,50 kDa

研究背景 (Background)


Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

产品名称: Annexin A7 Rabbit Monoclonal Antibody
产品货号: AMRe84099




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Annexin A7 expression in (1) Jurkat cell lysate; (2) Rat kidney lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:1K dilution)

注意事项 (Note)

For research use only .