

产品名称: HSPC150 Rabbit Monoclonal Antibody
产品货号: AMRe84108



产品概述 (Summary)

产品名称 (Production Name)	HSPC150 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPC150
别名 (Alternative Names)	HSPC150; PIG50; Ube2t; UBE2T
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NPD8.A synthesized peptide derived from human UBE2T

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	23 kDa

研究背景 (Background)

Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. Catalyzes


产品名称: HSPC150 Rabbit Monoclonal Antibody
产品货号: AMRe84108




monoubiquitination. Involved in mitomycin-C (MMC)-induced DNA repair. Acts as a specific E2 ubiquitin-conjugating enzyme for the Fanconi anemia complex by associating with E3 ubiquitin-protein ligase FANCL and catalyzing monoubiquitination of FANCD2, a key step in the DNA damage pathway.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of HSPC150 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .