

产品名称: EXOC2 Rabbit Monoclonal Antibody
产品货号: AMRe84109



产品概述 (Summary)

产品名称 (Production Name)	EXOC2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EXOC2
别名 (Alternative Names)	EXOC2; SEC5; SEC5L1; Sec5p;;Exocyst complex component 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96KP1.A synthesized peptide derived from human Exocyst complex component 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	104 kDa

研究背景 (Background)


产品名称: EXOC2 Rabbit Monoclonal Antibody
产品货号: AMRe84109




Component of the exocyst complex involved in the docking of exocytic vesicles with fusion sites on the plasma membrane.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of EXOC2 expression in BxPC-3 cell lysate.

产品名称: EXOC2 Rabbit Monoclonal Antibody
产品货号: AMRe84109



注意事项 (Note)

For research use only .