

产品名称: RAB22A Rabbit Monoclonal Antibody
产品货号: AMRe84110



产品概述 (Summary)

产品名称 (Production Name)	RAB22A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB22A
别名 (Alternative Names)	Rab22; RAB22A;; Rab22A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UL26.A synthesized peptide derived from human Rab22A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	22 kDa

研究背景 (Background)

Plays a role in endocytosis and intracellular protein transport. Mediates trafficking of TF from early endosomes to recycling


产品名称: RAB22A Rabbit Monoclonal Antibody
产品货号: AMRe84110




endosomes. Required for NGF-mediated endocytosis of NTRK1, and subsequent neurite outgrowth. Binds GTP and GDP and has low GTPase activity. Alternates between a GTP-bound active form and a GDP-bound inactive form.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




Western blot analysis of RAB22A expression in (1) MCF-7 cell lysate; (2) BxPC-3 cell lysate.


产品名称: RAB22A Rabbit Monoclonal Antibody
产品货号: AMRe84110




Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:1K dilution)



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.



产品名称: RAB22A Rabbit Monoclonal Antibody
产品货号: AMRe84110



Immunofluorescent analysis using the Antibody at 1:150 dilution.

注意事项 (Note)

For research use only .