

产品名称: PTBP2 Rabbit Monoclonal Antibody
产品货号: AMRe84113



产品概述 (Summary)

产品名称 (Production Name)	PTBP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PTBP2
别名 (Alternative Names)	MIBP; nPTB; nPTB5; nPTB6; nPTB7; nPTB8; PTB; PTBLP; Ptbp2;;PTBP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UKA9.A synthesized peptide derived from human PTBP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa ; Observed MW: 55 kDa

研究背景 (Background)

RNA-binding protein which binds to intronic polypyrimidine tracts and mediates negative regulation of exons splicing. May


产品名称: PTBP2 Rabbit Monoclonal Antibody
产品货号: AMRe84113




antagonize in a tissue-specific manner the ability of NOVA1 to activate exon selection. In addition to its function in pre-mRNA splicing, plays also a role in the regulation of translation.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of PTBP2 expression in Neuro2a cell lysate.

注意事项 (Note)

For research use only .