

产品名称: Palladin Rabbit Monoclonal Antibody
产品货号: AMRe84114



产品概述 (Summary)

产品名称 (Production Name)	Palladin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Palladin
别名 (Alternative Names)	CGI15; MYN; Palld; PNCA1; SIH002; Stx13bp1;;Palladin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8WX93.A synthesized peptide derived from human Palladin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 151 kDa ; Observed MW: 130 kDa

研究背景 (Background)

Cytoskeletal protein required for organization of normal actin cytoskeleton. Roles in establishing cell morphology, motility,


产品名称: Palladin Rabbit Monoclonal Antibody
产品货号: AMRe84114




cell adhesion and cell-extracellular matrix interactions in a variety of cell types. May function as a scaffolding molecule with the potential to influence both actin polymerization and the assembly of existing actin filaments into higher-order arrays. Binds to proteins that bind to either monomeric or filamentous actin. Localizes at sites where active actin remodeling takes place, such as lamellipodia and membrane ruffles.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



产品名称: Palladin Rabbit Monoclonal Antibody
产品货号: AMRe84114



Western blot analysis of Palladin expression in PC3 cell lysate.

注意事项 (Note)

For research use only .